



WASP[®]

EMPOWERED BY BIOMÉRIEUX

We adapt to you. You don't have to adapt to us.



WASP




PIONEERING DIAGNOSTICS



YOUR LAB. YOUR WAY. YOUR WASP®.

Build your workstation to suit your lab! WASP® is a **unique open-platform, modular system** that grows with you as your needs change. Combine it with services available only through **WASP® empowered by bioMérieux for ultimate automation.**

THE AGILE & FLEXIBLE AUTOMATED SPECIMEN PROCESSING SYSTEM



QUALITY FROM START TO FINISH

- As microbiologists, we understand the need to pamper specimens to get the best results. WASP® offers easy solutions for quality colony growth:
- Use with compatible **Liquid Based Media (LBM®)** and **Liquid Amies Collection SWAB (ESWAB®)** ranges for **specimen protection** throughout transport, **optimized sample conservation & release** and **streamlined workflow** with standardized sample containers
 - Fully validated compatibility with **bioMérieux's broad range of high quality pre-plated media (PPM)**

CUSTOMIZATION FOR OPTIMAL COLONY ISOLATION

- Keep the standardization of automation and team it with specimen inoculation flexibility to suit your lab needs.
- Total versatility with an extensive and customizable library of **classic streaking patterns** that reproduce human gestures
 - Choose your specimen transfer volumes: Loops for **1 µl, 10 µl or 30 µl** and even twin loops for efficient biplate streaking
 - Optional **loops sterilization between quadrants** for optimal colony isolation of samples with high bacterial loads
 - Work with **any specimen type** (no pipetting required)


THE MOST AUTOMATED SPECIMEN PROCESSOR AVAILABLE

WASP® fits into **your lab** and lets you pick and mix the automation add-ons you need – we adapt to you, you don't have to adapt to us!

Automated Gram slide preparation



Pipetor module from 10 µl to 900 µl



Automated broth preparation




WASP® CAN OPERATE 24/7



Continuous Loading



Urine Turntable



Antibiotic Disc Deposit system



Why choose WASP® empowered by bioMérieux ?

WASP® empowered by bioMérieux is, simply, more. We back your WASP® instrument with **flexible, innovative services** developed using our 50 years of microbiology expertise.


CONCRETELY, WE:

- Partner with you for the life of your WASP through our **5-star support service package**
- Provide relevant **training and education** by microbiology experts
- Leverage our **microbiology know-how** to optimize your WASP and its use
- **Smoothly integrate automation and informatics** facilitating the continuous improvement of your processes





WASP® FOR YOUR NEEDS TODAY.



WASP® IN THE WORLD

THE LARGEST STREAKER INSTALLED BASE


Please contact your local bioMérieux representative for availability of these solutions in your country.

Manufactured by
Via F. Perotti 10 • 25125 Brescia • Italy



Distributed by
bioMérieux S.A. • 69280 Marcy l'Etoile • France • Tel.: + 33 (0)4 78 87 20 00 • Fax: +33 (0)4 78 87 20 90
www.biomerieux.com

WASPLab®, AVAILABLE FOR YOUR LAB TOMORROW



Your need for extra productivity may evolve over time. No problem! Count on us to find the right set-up to suit your new needs. Your WASP® can be a launch pad for increased lab efficiency by evolving to WASPLab® over time.
It's a smart investment in the age of the automated lab.